

Analyte Test Name	Specimen Type	Specimen Stability (Ambient, Refrig, and Frozen)
Ethanol (Alcohol)	Urine	
Amphetamine	Urine	Ambient: 14 days, Refrigerated: 14 days
Methamphetamine	Urine	Ambient: 14 days, Refrigerated: 14 days
MDA	Urine	Ambient: 14 days, Refrigerated: 14 days
MDMA	Urine	Ambient: 14 days, Refrigerated: 14 days
MDEA	Urine	Ambient: 14 days, Refrigerated: 14 days
Amobarbital/Pentobarbital	Urine	Ambient: 12 days, Refrigerated: 12 days
Butalbital	Urine	Ambient: 12 days, Refrigerated: 12 days
Phenobarbital	Urine	Ambient: 12 days, Refrigerated: 12 days
Secobarbital	Urine	Ambient: 12 days, Refrigerated: 12 days
Lorazepam	Urine	Ambient: 7 days, Refrigerated: 14 days

Oxazepam	Urine	Ambient: 14 days, Refrigerated: 14 days
Nordiazepam	Urine	Ambient: 14 days, Refrigerated: 14 days
Temazepam	Urine	Ambient: 14 days, Refrigerated: 14 days
Nitrazepam Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Clonazepam Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Flunitrazepam Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Triazolam Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Flurazepam Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Alprazolam Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Midazolam Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Estazolam	Urine	Ambient: 14 days, Refrigerated: 14 days

Marijuana Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Cocaine Metabolites	Urine	Ambient: 14 days, Refrigerated: 14 days
Nicotine Metabolite	Urine Urine (Clinical)	Ambient: 12 days, Refrigerated: 12 days
Methadone	Urine	Ambient: 14 days, Refrigerated: 14 days
Methadone Metabolite (EDDP)	Urine	Ambient: 14 days, Refrigerated: 14 days
Phencyclidine (PCP)	Urine	Ambient: 14 days, Refrigerated: 14 days
ADB-BUTINACA N-(4-hydroxybutyl) mtb	Urine	Ambient: 14 days Refrigerated: 14 days
AMP-4en-PINACA	Urine	Ambient: 14 days, Refrigerated: 14 days
4F MDMB-BUTINACA N-butanoic acid mtb	Urine	Ambient: 14 days, Refrigerated: 14 days
4F MDMB-BUTICA butanoic acid mtb	Urine	Ambient: 14 days, Refrigerated: 14 days
ADB-BINACA	Urine	Ambient: 14 days, Refrigerated: 14 days

ADB-BUTINACA	Urine	Ambient: 14 days, Refrigerated: 14 days
5F MDMB-PICA mtb 7	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 018 N-pentanoic acid mtb	Urine	Ambient: 14 days, Refrigerated: 14 days
ADB-4en-PINACA	Urine	Ambient: 14 days, Refrigerated: 14 days
5F ADB mtb 7	Urine	Ambient: 14 days, Refrigerated: 14 days
4-cyano CUMYL-BUTINACA	Urine	Ambient: 14 days, Refrigerated: 14 days
4F MDMB-BUTICA	Urine	Ambient: 14 days, Refrigerated: 14 days
4F ABUTINACA N-(4-hydroxybutyl) mtb	Urine	Ambient: 14 days, Refrigerated: 14 days
MDMB-4en-PINACA butanoic acid mtb	Urine	Ambient: 14 days, Refrigerated: 14 days
4F MDMB-BUTINACA	Urine	Ambient: 14 days, Refrigerated: 14 days
5F MDMB-PICA	Urine	Ambient: 14 days, Refrigerated: 14 days

ADB-HEXINACA	Urine	Ambient: 14 days, Refrigerated: 14 days
MDMB-CHMICA mtb M2	Urine	Ambient: 14 days, Refrigerated: 14 days
Ketamine Metabolite	Urine	Ambient:14 days, Refrigerated: 14 days
Zolpidem	Urine	Ambient: 14 days, Refrigerated: 14 days
Zolpidem Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Gabapentin	Urine	Ambient: 14 days, Refrigerated: 14 days
Pregabalin	Urine	Ambient: 14 days, Refrigerated: 14 days
Tapentadol	Urine	Ambient: 14 days, Refrigerated: 14 days
Diphenhydramine	Urine	Ambient: 12 days, Refrigerated: 12 days
Chlorpheniramine	Urine	Ambient: 12 days, Refrigerated: 12 days
Doxylamine	Urine	Ambient: 12 days, Refrigerated: 12 days

Pheniramine	Urine	Ambient: 12 days, Refrigerated: 12 days
-------------	-------	--

Brompheniramine	Urine	Ambient: 12 days, Refrigerated: 12 days
-----------------	-------	--

Pseudo/Ephedrine	Urine	Ambient: 12 days, Refrigerated: 12 days
------------------	-------	--

Phentermine	Urine	Ambient: 12 days, Refrigerated: 12 days
-------------	-------	--

Phenylpropanolamine	Urine	Ambient: 12 days, Refrigerated: 12 days
---------------------	-------	--

Diethylpropion	Urine	Ambient: 12 days, Refrigerated: 12 days
----------------	-------	--

Phenmetrazine	Urine	Ambient: 12 days, Refrigerated: 12 days
---------------	-------	--

Phendimetrazine	Urine	Ambient: 12 days, Refrigerated: 12 days
-----------------	-------	--

Naloxone	Urine	Ambient: 14 days, Refrigerated: 14 days
----------	-------	--

Norbuprenorphine	Urine	Ambient: 14 days, Refrigerated: 14 days
------------------	-------	--

Buprenorphine	Urine	Ambient: 14 days, Refrigerated: 14 days
---------------	-------	--

Eszopiclone/Zopiclone	Urine	Ambient: 12 days, Refrigerated: 12 days
Morphine	Urine	Ambient: 14 days, Refrigerated: 14 days
Hydromorphone	Urine	Ambient: 14 days, Refrigerated: 14 days
Dihydrocodeine	Urine	Ambient: 14 days, Refrigerated: 14 days
Codeine	Urine	Ambient: 14 days, Refrigerated: 14 days
Hydrocodone	Urine	Ambient: 14 days, Refrigerated: 14 days
Norhydrocodone	Urine	Ambient: 14 days, Refrigerated: 14 days
6-Acetylmorphine (Heroin)	Urine	Ambient: 14 days, Refrigerated: 14 days
Propofol Glucuronide	Urine	Ambient: >14 days, Refrigerated: >14 days
Buspirone	Urine	Ambient: 14 days, Refrigerated: 14 days
Buspirone Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days

Clomipramine

Urine

Ambient: 14 days,
Refrigerated: 14 days

Clomipramine Metabolite

Urine

Ambient: 14 days,
Refrigerated: 14 days

Lamotrigine

Urine

Ambient: 14 days,
Refrigerated: 14 days

Acetyl Fentanyl

Urine

Ambient: 14 days,
Refrigerated: 14 days

Mephedrone

Urine

Ambient: 14 days,
Refrigerated: 14 days

MDPV

Urine

Ambient: 14 days,
Refrigerated: 14 days

Methylone

Urine

Ambient: 14 days,
Refrigerated: 14 days

Methedrone

Urine

Ambient: 14 days,
Refrigerated: 14 days

Pentylone

Urine

Ambient: 14 days,
Refrigerated: 14 days

Alpha-PVP

Urine

Ambient: 14 days,
Refrigerated: 14 days

Pyrovalerone

Urine

Ambient: 14 days,
Refrigerated: 14 days

2C-E(Europa)	Urine	Ambient: 12 days, Refrigerated: 12 days
Mitragynine (Kratom)	Urine	Ambient: 14 days, Refrigerated: 14 days
Norfluoxetine	Urine	Ambient: 14 days, Refrigerated: 14 days
Sertraline	Urine	Ambient: 14 days, Refrigerated: 14 days
Duloxetine	Urine	Ambient: 14 days, Refrigerated: 14 days
Fluoxetine	Urine	Ambient: 14 days, Refrigerated: 14 days
Nortriptyline	Urine	Ambient: 14 days, Refrigerated: 14 days
Amitriptyline	Urine	Ambient: 14 days, Refrigerated: 14 days
Nordoxepin	Urine	Ambient: 14 days, Refrigerated: 14 days
Desipramine	Urine	Ambient: 14 days, Refrigerated: 14 days
Methylphenidate Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days

Oxymorphone

Urine

Ambient: 14 days,
Refrigerated: 14 days

Oxycodone

Urine

Ambient: 14 days,
Refrigerated: 14 days

Noroxycodone

Urine

Ambient: 14 days,
Refrigerated: 14 days

Tramadol

Urine

Ambient: 14 days,
Refrigerated: 14 days

O-Desmethyl Tramadol

Urine

Ambient: 14 days,
Refrigerated: 14 days

Meprobamate

Urine

Ambient: 14 days,
Refrigerated: 14 days

Meperidine Metabolite

Urine

Ambient: 14 days,
Refrigerated: 14 days

Fentanyl Metabolite

Urine

Ambient: 14 days,
Refrigerated: 14 days

Acetyl Fentanyl

Urine

Ambient: 14 days,
Refrigerated: 14 days

Fentanyl

Urine

Ambient: 14 days,
Refrigerated: 14 days

Suvorexant

Urine

Ambient: 14 days,
Refrigerated: 14 days

Zaleplon	Urine	Ambient: 14 days, Refrigerated: 14 days
Baclofen	Urine	Ambient: 14 days, Refrigerated: 14 days
Cyclobenzaprine	Urine	Ambient: 14 days, Refrigerated: 14 days
Methocarbamol	Urine	Ambient: 14 days, Refrigerated: 14 days
Tizanidine	Urine	Ambient: 14 days, Refrigerated: 14 days
Aripiprazole	Urine	Ambient: 14 days, Refrigerated: 14 days
Aripiprazole/Brexpiprazole Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Chlorpromazine	Urine	Ambient: 2-3 days, Refrigerated: 14 days
Flupenthixol	Urine	Ambient: 2-3 days, Refrigerated: 14 days
Clozapine	Urine	Ambient: 14 days, Refrigerated: 14 days
Clozapine Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days

Haloperidol	Urine	Ambient: 14 days, Refrigerated: 14 days
Lurasidone	Urine	Ambient: 2-3 days, Refrigerated: 14 days
Molindone	Urine	Ambient: 14 days, Refrigerated: 14 days
Paliperidone	Urine	Ambient: 14 days, Refrigerated: 14 days
Perphenazine	Urine	Ambient: Not Stable, Refrigerated: 2-3 days
Quetiapine	Urine	Ambient: 14 days, Refrigerated: 14 days
Risperidone	Urine	Ambient: 14 days, Refrigerated: 14 days
Thioridazine	Urine	Ambient: Not Stable, Refrigerated: 2-3 days
Thiothixene	Urine	Ambient: 7 days, Refrigerated: 14 days
Trifluoperazine	Urine	Ambient: 2-3 days, Refrigerated: 2-3 days
Amoxapine	Urine	Ambient: 14 days, Refrigerated: 14 days

Amoxapine Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Clomipramine	Urine	Ambient: 14 days, Refrigerated: 14 days
Clomipramine Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Mirtazapine	Urine	Ambient: 14 days, Refrigerated: 14 days
Trimipramine	Urine	Ambient: 14 days, Refrigerated: 14 days
Tranlycypromine	Urine	Ambient: 14 days, Refrigerated: 14 days
Bupropion	Urine	Ambient: 2-3 days, Refrigerated: 14 days
Bupropion Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Trazodone	Urine	Ambient: 14 days, Refrigerated: 14 days
Venlafaxine Metabolite	Urine	Ambient: 14 days, Refrigerated: 14 days
Venlafaxine	Urine	Ambient: 14 days, Refrigerated: 14 days

Paroxetine	Urine	Ambient: 14 days, Refrigerated: 14 days
------------	-------	--

Citalopram/Escitalopram	Urine	Ambient: 14 days, Refrigerated: 14 days
-------------------------	-------	--

Levorphanol/Dextrophan	Urine	Ambient: 14 days, Refrigerated: 14 days
------------------------	-------	--

Dextromethorphan	Urine	Ambient: 14 days, Refrigerated: 14 days
------------------	-------	--

Naltrexone	Urine	Ambient: 14 days, Refrigerated: 14 days
------------	-------	--

Nalbuphine	Urine	Ambient: 14 days, Refrigerated: 14 days
------------	-------	--

Pentazocine	Urine	Ambient: 14 days, Refrigerated: 14 days
-------------	-------	--

Butorphanol	Urine	Ambient: 14 days, Refrigerated: 14 days
-------------	-------	--

Norpropoxyphene	Urine	Ambient: 12 days, Refrigerated: 12 days
-----------------	-------	--

JWH 073 hydroxybutyl	Urine	Ambient: 14 days, Refrigerated: 14 days
----------------------	-------	--

JWH 073 butanoic acid	Urine	Ambient: 14 days, Refrigerated: 14 days
-----------------------	-------	--

JWH 018 pentanoic acid	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 018 hydroxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 073 hydroxyindol	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 018 hydroxyindol	Urine	Ambient: 14 days, Refrigerated: 14 days
RCS4 N-5-carboxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 250-N-5-carboxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
AM2201 N-4-hydroxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 081-N-5-hydroxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 019-N-6-hydroxyhexyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 122 pentanoic acid	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 210-N-5-carboxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days

XLR11 hydroxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
UR144 pentanoic acid	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 073 hydroxybutyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 073 butanoic acid	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 018 pentanoic acid	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 018 hydroxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 073 hydroxyindol	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 018 hydroxyindol	Urine	Ambient: 14 days, Refrigerated: 14 days
RCS4 N-5-carboxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 250-N-5-carboxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
AM2201 N-4-hydroxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days

JWH 081-N-5-hydroxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 019-N-6-hydroxyhexyl	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 122 pentanoic acid	Urine	Ambient: 14 days, Refrigerated: 14 days
JWH 210-N-5-carboxypentyl	Urine	Ambient: 14 days, Refrigerated: 14 days
XLR11 hydroxypentyl	Urine	Ambient: 7 days, Refrigerated: 14 days
UR144 pentanoic acid	Urine	Ambient: 7 days, Refrigerated: 14 days
LSD	Urine	Ambient: 14 days, Refrigerated: 14 days
OxoHydroxyLSD	Urine	Ambient: 1 days, Refrigerated: 2 days
Norbuprenorphine	Urine	Ambient: 14 days, Refrigerated: 14 days
Buprenorphine	Urine	Ambient: 14 days, Refrigerated: 14 days
Ethylsulfate	Urine	Ambient: 7 days, Refrigerated: 7 days

Ethylglucuronide	Urine	Ambient: 7 days, Refrigerated: 7 days
Nicotine Metabolite	Urine (Clinical	Ambient: 12 days, Refrigerated: 12 days
Amphetamine	Oral Fluid	Ambient: 14 days
Methamphetamine	Oral Fluid	Ambient: 14 days
Amobarbital	Oral Fluid	Ambient: 14 days, Refrigerated: 14 days
Butalbital	Oral Fluid	Ambient: 14 days, Refrigerated: 14 days
Pentobarbital	Oral Fluid	Ambient: 14 days, Refrigerated: 14 days
Phenobarbital	Oral Fluid	Ambient: 14 days, Refrigerated: 14 days
Secobarbital	Oral Fluid	Ambient: 14 days, Refrigerated: 14 days
Lorazepam	Oral Fluid	Ambient: 14 days
Oxazepam	Oral Fluid	Ambient: 14 days
Clonazepam	Oral Fluid	Ambient: 14 days

Bromazepam

Oral Fluid

Ambient: 14 days

Nitrazepam

Oral Fluid

Ambient: 14 days

Midazolam

Oral Fluid

Ambient: 14 days

Nordiazepam

Oral Fluid

Ambient: 14 days

Flunitrazepam

Oral Fluid

Ambient: < 7 days

Lormetazepam

Oral Fluid

Ambient: 14 days

Triazolam

Oral Fluid

Ambient: 14 days

Estazolam

Oral Fluid

Ambient: 14 days

Temazepam

Oral Fluid

Ambient: 14 days

Alprazolam

Oral Fluid

Ambient: 14 days

Diazepam

Oral Fluid

Ambient: 14 days

Prazepam	Oral Fluid	Ambient: 14 days
----------	------------	------------------

Clonazepam Metabolite	Oral Fluid	Ambient: 14 days
-----------------------	------------	------------------

Cocaine	Oral Fluid	Ambient: 7-14 days
---------	------------	--------------------

Benzoyllecgonine	Oral Fluid	Ambient: 14 days
------------------	------------	------------------

Buprenorphine	Oral Fluid	Ambient: 14 days
---------------	------------	------------------

Norbuprenorphine	Oral Fluid	Ambient: 14 days
------------------	------------	------------------

Fentanyl Metabolite	Oral Fluid	Ambient: 14 days
---------------------	------------	------------------

Fentanyl	Oral Fluid	Ambient: 14 days
----------	------------	------------------

Meperidine Metabolite	Oral Fluid	Ambient: 14 days
-----------------------	------------	------------------

Meperidine	Oral Fluid	Ambient: 14 days
------------	------------	------------------

Tramadol	Oral Fluid	Ambient: 14 days
----------	------------	------------------

Meprobamate	Oral Fluid	Ambient: 14 days
-------------	------------	------------------

Carisoprodol	Oral Fluid	Ambient: 14 days
--------------	------------	------------------

Methadone	Oral Fluid	Ambient: 7-14 days
-----------	------------	--------------------

Methadone Metabolite (EDDP)	Oral Fluid	Ambient: 14 days
-----------------------------	------------	------------------

Morphine	Oral Fluid	Ambient: 14 days
----------	------------	------------------

Hydromorphone	Oral Fluid	Ambient: 14 days
---------------	------------	------------------

Dihydrocodeine	Oral Fluid	Ambient: 14 days
----------------	------------	------------------

Codeine	Oral Fluid	Ambient: 14 days
---------	------------	------------------

6-Acetylmorphine (Heroin)	Oral Fluid	Ambient: 14 days
---------------------------	------------	------------------

Hydrocodone	Oral Fluid	Ambient: 14 days
-------------	------------	------------------

Oxymorphone	Oral Fluid	Ambient: 14 days
-------------	------------	------------------

Oxycodone

Oral Fluid

Ambient: 14 days

Phencyclidine (PCP)

Oral Fluid

Ambient: 14 days

Delta-9-THC

Oral Fluid

N/A

Carboxy-THC (metabolite)

Oral Fluid

Ambient: 14 days,
Refrigerated: 14 days